




HA Antibody

Product Code	CSB-PA607773LA01IER
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0HD60
Immunogen	Recombinant Influenza A virus Hemagglutinin protein (18-343AA)
Raised In	Rabbit
Species Reactivity	Influenza A virus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemagglutinin [Cleaved into: Hemagglutinin HA1 chain; Hemagglutinin HA2 chain], HA
Immunogen Species	Influenza A virus
Research Area	Others
Target Names	HA

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: HA antibody at 3μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 64 kDa

observed band size: 64 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.