



PRDM2 Antibody, Biotin conjugated

Product Code	CSB-PA613379LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13029
Immunogen	Recombinant Human PR domain zinc finger protein 2 protein (733-961AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PR domain zinc finger protein 2 (EC 2.1.1.43) (GATA-3-binding protein G3B) (Lysine N-methyltransferase 8) (MTB-ZF) (MTE-binding protein) (PR domain-containing protein 2) (Retinoblastoma protein-interacting zinc finger protein) (Zinc finger protein RIZ), PRDM2, KMT8 RIZ
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PRDM2