





DYNC1H1 Antibody

Product Code	CSB-PA613491LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14204
Immunogen	Recombinant Human Cytoplasmic dynein 1 heavy chain 1 protein (12-66AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytoplasmic dynein 1 heavy chain 1 (Cytoplasmic dynein heavy chain 1) (Dynein heavy chain, cytosolic), DYNC1H1, DHC1 DNCH1 DNCL DNECL DYHC KIAA0325
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DYNC1H1

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA613491LA01HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA613491LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.