



DYNC1H1 Antibody, Biotin conjugated

Product Code	CSB-PA613491LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14204
Immunogen	Recombinant Human Cytoplasmic dynein 1 heavy chain 1 protein (12-66AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytoplasmic dynein 1 heavy chain 1 (Cytoplasmic dynein heavy chain 1) (Dynein heavy chain, cytosolic), DYNC1H1, DHC1 DNCH1 DNCL DNECL DYHC KIAA0325
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DYNC1H1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.