



KIR2DS3 Antibody, HRP conjugated

Product Code	CSB-PA613530LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14952
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 2DS3 protein (102-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Killer cell immunoglobulin-like receptor 2DS3 (MHC class I NK cell receptor) (Natural killer-associated transcript 7) (NKAT-7), KIR2DS3, NKAT7
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KIR2DS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.