



KIR2DS3 Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA613530LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q14952 |
| Immunogen | Recombinant Human Killer cell immunoglobulin-like receptor 2DS3 protein (102-241AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Killer cell immunoglobulin-like receptor 2DS3 (MHC class I NK cell receptor) (Natural killer-associated transcript 7) (NKAT-7), KIR2DS3, NKAT7 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | KIR2DS3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |