



SLAMF1 Antibody, FITC conjugated

Product Code	CSB-PA614403LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13291
Immunogen	Recombinant Human Signaling lymphocytic activation molecule protein (69-204AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signaling lymphocytic activation molecule (CDw150) (IPO-3) (SLAM family member 1) (CD antigen CD150), SLAMF1, SLAM
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SLAMF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.