



PLD1 Antibody, FITC conjugated

Product Code	CSB-PA614406LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13393
Immunogen	Recombinant Human Phospholipase D1 protein (485-627AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phospholipase D1 (PLD 1) (hPLD1) (EC 3.1.4.4) (Choline phosphatase 1) (Phosphatidylcholine-hydrolyzing phospholipase D1), PLD1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PLD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.