



PTK6 Antibody, Biotin conjugated

Product Code	CSB-PA614448LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13882
Immunogen	Recombinant Human Protein-tyrosine kinase 6 protein (1-113AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-tyrosine kinase 6 (EC 2.7.10.2) (Breast tumor kinase) (Tyrosine-protein kinase BRK), PTK6, BRK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTK6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.