



CAMK1 Antibody, FITC conjugated

Product Code	CSB-PA614509LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14012
Immunogen	Recombinant Human Calcium/calmodulin-dependent protein kinase type 1 protein (166-265AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium/calmodulin-dependent protein kinase type 1 (EC 2.7.11.17) (CaM kinase I) (CaM-KI) (CaM kinase I alpha) (CaMKI-alpha), CAMK1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CAMK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.