



SSX1 Antibody, HRP conjugated

Product Code	CSB-PA614964LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16384
Immunogen	Recombinant Human Protein SSX1 protein (112-188AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein SSX1 (Cancer/testis antigen 5.1) (CT5.1) (Synovial sarcoma, X breakpoint 1), SSX1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SSX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.