



MAP2K5 Antibody, FITC conjugated

Product Code	CSB-PA615667HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13163
Immunogen	Recombinant Human Dual specificity mitogen-activated protein kinase kinase 5 protein (1-190AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity mitogen-activated protein kinase kinase 5 (MAP kinase kinase 5) (MAPKK 5) (EC 2.7.12.2) (MAPK/ERK kinase 5) (MEK 5), MAP2K5, MEK5 MKK5 PRKMK5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP2K5