



KLRD1 Antibody, FITC conjugated

Product Code	CSB-PA615672LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13241
Immunogen	Recombinant Human Natural killer cells antigen CD94 protein (42-168AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Natural killer cells antigen CD94 (KP43) (Killer cell lectin-like receptor subfamily D member 1) (NK cell receptor) (CD antigen CD94), KLRD1, CD94
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KLRD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.