



KLRD1 Antibody, Biotin conjugated

Product Code	CSB-PA615672LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13241
Immunogen	Recombinant Human Natural killer cells antigen CD94 protein (42-168AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Natural killer cells antigen CD94 (KP43) (Killer cell lectin-like receptor subfamily D member 1) (NK cell receptor) (CD antigen CD94), KLRD1, CD94
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KLRD1