



UAP1 Antibody, Biotin conjugated

Product Code	CSB-PA618067LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16222
Immunogen	Recombinant Human UDP-N-acetylhexosamine pyrophosphorylase protein (1-522AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-N-acetylhexosamine pyrophosphorylase (Antigen X) (AGX) (Sperm-associated antigen 2) [Includes: UDP-N-acetylgalactosamine pyrophosphorylase (EC 2.7.7.83) (AGX-1); UDP-N-acetylglucosamine pyrophosphorylase (EC 2.7.7.23) (AGX-2)], UAP1, SPAG2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	UAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.