



UGT8 Antibody

Product Code	CSB-PA618101LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16880
Immunogen	Recombinant Human 2-hydroxyacylsphingosine 1-beta-galactosyltransferase protein (38-402AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	2-hydroxyacylsphingosine 1-beta-galactosyltransferase (EC 2.4.1.47) (Ceramide UDP-galactosyltransferase) (Cerebroside synthase) (UDP-galactose-ceramide galactosyltransferase), UGT8, CGT UGT4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	UGT8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.