



DGKZ Antibody, Biotin conjugated

Product Code	CSB-PA618788LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13574
Immunogen	Recombinant Human Diacylglycerol kinase zeta protein (11-152AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diacylglycerol kinase zeta (DAG kinase zeta) (EC 2.7.1.107) (Diglyceride kinase zeta) (DGK-zeta), DGKZ, DAGK6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DGKZ