



SCARB2 Antibody, FITC conjugated

Product Code	CSB-PA619859LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14108
Immunogen	Recombinant Human Lysosome membrane protein 2 protein (47-128AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysosome membrane protein 2 (85 kDa lysosomal membrane sialoglycoprotein) (LGP85) (CD36 antigen-like 2) (Lysosome membrane protein II) (LIMP II) (Scavenger receptor class B member 2) (CD antigen CD36), SCARB2, CD36L2 LIMP2 LIMPII
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SCARB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.