



UBE3C Antibody, Biotin conjugated

Product Code	CSB-PA621871LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15386
Immunogen	Recombinant Human Ubiquitin-protein ligase E3C protein (310-377AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-protein ligase E3C (EC 2.3.2.26) (HECT-type ubiquitin transferase E3C) (HectH2), UBE3C, KIAA0010 KIAA10
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBE3C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.