



TRADD Antibody, HRP conjugated

Product Code	CSB-PA621879EB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15628
Immunogen	Recombinant Human Tumor necrosis factor receptor type 1-associated DEATH domain protein (1-312AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor type 1-associated DEATH domain protein (TNFR1-associated DEATH domain protein) (TNFRSF1A-associated via death domain), TRADD
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRADD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.