



# FOXF2 Antibody

Product Code	CSB-PA623661LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12947
Immunogen	Recombinant Human Forkhead box protein F2 protein (226-364AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein F2 (Forkhead-related activator 2) (FREAC-2) (Forkhead-related protein FKHL6) (Forkhead-related transcription factor 2), FOXF2, FKHL6 FREAC2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FOXF2