





MDC1 Antibody

Product Code	CSB-PA623929LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14676
Immunogen	Recombinant Human Mediator of DNA damage checkpoint protein 1 protein (1892-2082AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mediator of DNA damage checkpoint protein 1 (Nuclear factor with BRCT domains 1), MDC1, KIAA0170 NFBD1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MDC1

Image



Immunohistochemistry analysis of human breast cancer using CSB-PA623929LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA623929LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.