



MDC1 Antibody, Biotin conjugated

Product Code	CSB-PA623929LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14676
Immunogen	Recombinant Human Mediator of DNA damage checkpoint protein 1 protein (1892-2082AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mediator of DNA damage checkpoint protein 1 (Nuclear factor with BRCT domains 1), MDC1, KIAA0170 NFBD1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MDC1