



vjbR Antibody, HRP conjugated

Product Code	CSB-PA651339YB01BMN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2YJ50
Immunogen	Recombinant Brucella abortus HTH-type quorum sensing-dependent transcriptional regulator VjbR protein (1-259AA)
Raised In	Rabbit
Species Reactivity	Brucella abortus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HTH-type quorum sensing-dependent transcriptional regulator VjbR, vjbR
Immunogen Species	Brucella abortus (strain 2308)
Research Area	Others
Target Names	vjbR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.