



vacA Antibody, HRP conjugated

Product Code	CSB-PA669657LB01HUS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q48253
Immunogen	Recombinant Helicobacter pylori Vacuolating cytotoxin autotransporter protein (1038-1310AA)
Raised In	Rabbit
Species Reactivity	Helicobacter pylori
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolating cytotoxin autotransporter [Cleaved into: Vacuolating cytotoxin; Vacuolating cytotoxin translocator], vacA
Immunogen Species	Helicobacter pylori (Campylobacter pylori)
Research Area	Others
Target Names	vacA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.