



LEXM Antibody, FITC conjugated

Product Code	CSB-PA670268LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3ZCV2
Immunogen	Recombinant Human Lymphocyte expansion molecule protein (1-418AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphocyte expansion molecule, LEXM, LEM
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LEXM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.