




Pollen allergen Phl p 5b Antibody

Product Code	CSB-PA671435HA01EUQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q40963
Immunogen	Recombinant Phleum pratense Pollen allergen Phl p 5b protein (20-284AA)
Raised In	Rabbit
Species Reactivity	Phleum pRatense
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pollen allergen Phl p 5b (Allergen Phl p Vb) (allergen Phl p 5b) (Fragment)
Immunogen Species	Phleum pratense
Research Area	Others

Image



Western Blot

Positive WB detected in Recombinant protein
All lanes: Pollen allergen Phl p 5b antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
predicted band size: 43 kDa
observed band size: 43 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.