



speB Antibody, HRP conjugated

Product Code	CSB-PA676488LB01SBAF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q48R29
Immunogen	Recombinant Streptococcus pyogenes serotype M28 Streptopain protein (146-398AA)
Raised In	Rabbit
Species Reactivity	Streptococcus pyogenes serotype M28
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Streptopain (EC 3.4.22.10) (Exotoxin type B) (SPE B) (Streptococcal cysteine proteinase) (Streptococcus peptidase A) (SPP), speB
Immunogen Species	Streptococcus pyogenes serotype M28 (strain MGAS6180)
Research Area	Others
Target Names	speB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.