



pat Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA701139HD01FOQ |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q57146 |
| Immunogen | Recombinant Streptomyces viridochromogenes Phosphinothricin N-acetyltransferase protein (1-183AA) |
| Raised In | Rabbit |
| Species Reactivity | Streptomyces viridochromogenes |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Phosphinothricin N-acetyltransferase (PPT N-acetyltransferase) (EC 2.3.1.183) (Phosphinothricin-resistance protein), pat |
| Immunogen Species | Streptomyces viridochromogenes (strain DSM 40736 / JCM 4977 / BCRC 1201 / Tue 494) |
| Research Area | Others |
| Target Names | pat |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |