




pat Antibody

Product Code	CSB-PA701139YA01FOQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q57146
Immunogen	Recombinant Streptomyces viridochromogenes Phosphinothricin N-acetyltransferase protein (1-183AA)
Raised In	Rabbit
Species Reactivity	Streptomyces viridochromogenes
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptomyces viridochromogenes
Research Area	Others
Target Names	pat

Image



Western Blot

Positive WB detected in: pat 293T Transfected lysate, 293T non-Transfected lysate

All lanes: pat antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 49 kDa

Observed band size: 49 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.