



pat Antibody, Biotin conjugated

Product Code	CSB-PA701139YD01FOQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q57146
Immunogen	Recombinant Streptomyces viridochromogenes Phosphinothricin N-acetyltransferase protein (1-183AA)
Raised In	Rabbit
Species Reactivity	Streptomyces viridochromogenes
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphinothricin N-acetyltransferase (PPT N-acetyltransferase) (EC 2.3.1.183) (Phosphinothricin-resistance protein), pat
Immunogen Species	Streptomyces viridochromogenes (strain DSM 40736 / JCM 4977 / BCRC 1201 / Tue 494)
Research Area	Others
Target Names	pat
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.