



nosZ Antibody

Product Code	CSB-PA703536LA01PBZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q51705
Immunogen	Recombinant Paracoccus denitrificans Nitrous-oxide reductase protein (58-553AA)
Raised In	Rabbit
Species Reactivity	Paracoccus denitrificans
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nitrous-oxide reductase (EC 1.7.2.4) (N(2)OR) (N2O reductase), nosZ
Immunogen Species	Paracoccus denitrificans
Research Area	Others
Target Names	nosZ