



nosZ Antibody, HRP conjugated

Product Code	CSB-PA703536LB01PBZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q51705
Immunogen	Recombinant Paracoccus denitrificans Nitrous-oxide reductase protein (58-553AA)
Raised In	Rabbit
Species Reactivity	Paracoccus denitrificans
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nitrous-oxide reductase (EC 1.7.2.4) (N(2)OR) (N2O reductase), nosZ
Immunogen Species	Paracoccus denitrificans
Research Area	Others
Target Names	nosZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.