



pvdA Antibody, HRP conjugated

Product Code	CSB-PA706553HB01EZ X
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q51548
Immunogen	Recombinant Pseudomonas aeruginosa L-ornithine N(5)-monooxygenase protein (1-443AA)
Raised In	Rabbit
Species Reactivity	Pseudomonas aeruginosa
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-ornithine N(5)-monooxygenase (EC 1.14.13.195) (L-ornithine N(5)-hydroxylase) (Ornithine hydroxylase) (L-ornithine N(5)-oxygenase) (Pyoverdin biosynthesis protein A), pvdA, pvd-1
Immunogen Species	Pseudomonas aeruginosa (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1)
Research Area	Others
Target Names	pvdA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.