




FHL5 Antibody

Product Code	CSB-PA725899LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TD97
Immunogen	Recombinant Human Four and a half LIM domains protein 5 protein (16-284AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Four and a half LIM domains protein 5 (FHL-5) (Activator of cAMP-responsive element modulator in testis) (Activator of CREM in testis), FHL5, ACT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FHL5

Image



Immunofluorescent analysis of HeLa cells using CSB-PA725899LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.