



SZT2 Antibody, FITC conjugated

Product Code	CSB-PA735962LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T011
Immunogen	Recombinant Human Protein SZT2 protein (2396-2517AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	KICSTOR complex protein SZT2 (Seizure threshold 2 protein homolog), SZT2, C1orf84 KIAA0467
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SZT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.