



hlgC Antibody, Biotin conjugated

Product Code	CSB-PA739984LD01SKV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6GE13
Immunogen	Recombinant Staphylococcus aureus Gamma-hemolysin component C protein (30-315AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-hemolysin component C, hlgC
Immunogen Species	Staphylococcus aureus (strain MRSA252)
Research Area	Others
Target Names	hlgC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.