



ATG14 Antibody, FITC conjugated

Product Code	CSB-PA741052LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZNE5
Immunogen	Recombinant Human Beclin 1-associated autophagy-related key regulator protein (1-492AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beclin 1-associated autophagy-related key regulator (Barkor) (Autophagy-related protein 14-like protein) (Atg14L), ATG14, ATG14L KIAA0831
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ATG14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.