



ATG14 Antibody, FITC conjugated

Product Code	CSB-PA741052LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZNE5
Immunogen	Recombinant Human Beclin 1-associated autophagy-related key regulator protein (1-492AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beclin 1-associated autophagy-related key regulator (Barkor) (Autophagy-related protein 14-like protein) (Atg14L), ATG14, ATG14L KIAA0831
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ATG14