



DNAH12 Antibody, Biotin conjugated

Product Code	CSB-PA747792LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZR08
Immunogen	Recombinant Human Dynein heavy chain 12, axonemal protein (2423-2632AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dynein heavy chain 12, axonemal (Axonemal beta dynein heavy chain 12) (Axonemal dynein heavy chain 12-like protein) (Axonemal dynein heavy chain 7-like protein) (Ciliary dynein heavy chain 12) (Dynein heavy chain 7-like, axonemal) (Dynein heavy chain domain-containing protein 2), DNAH12, DHC3 DLP12 DNAH12L DNAH7L DNAHC3 DNHD2 HDHC3 HL19
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNAH12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.