





NIPBL Antibody

Product Code	CSB-PA750379LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6KC79
Immunogen	Recombinant Human Nipped-B-like protein (2523-2697AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nipped-B-like protein (Delangin) (SCC2 homolog), NIPBL, IDN3 SCC2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NIPBL

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA750379LA01HU at dilution of 1:100



Western Blot

Positive WB detected in Recombinant protein

All lanes: NIPBL antibody at 4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 50 kDa

observed band size: 50 kDa



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.