



NIPBL Antibody, HRP conjugated

Product Code	CSB-PA750379LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6KC79
Immunogen	Recombinant Human Nipped-B-like protein (2523-2697AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nipped-B-like protein (Delangin) (SCC2 homolog), NIPBL, IDN3 SCC2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NIPBL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.