



TTC27 Antibody, HRP conjugated

Product Code	CSB-PA757569LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P3X3
Immunogen	Recombinant Human Tetratricopeptide repeat protein 27 protein (590-817AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 27 (TPR repeat protein 27), TTC27
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TTC27
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.