



# TTC27 Antibody, Biotin conjugated

Product Code	CSB-PA757569LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P3X3
Immunogen	Recombinant Human Tetratricopeptide repeat protein 27 protein (590-817AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 27 (TPR repeat protein 27), TTC27
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TTC27