



CNKSR3 Antibody, FITC conjugated

Product Code	CSB-PA757602LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P9H4
Immunogen	Recombinant Human Connector enhancer of kinase suppressor of ras 3 protein (301-555AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Connector enhancer of kinase suppressor of ras 3 (Connector enhancer of KSR 3) (CNK homolog protein 3) (CNK3) (CNKSR family member 3) (Maguin-like protein), CNKSR3, MAGI1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CNKSR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.