



CRTC1 Antibody, HRP conjugated

Product Code	CSB-PA757785LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UUV9
Immunogen	Recombinant Human CREB-regulated transcription coactivator 1 protein (1-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CREB-regulated transcription coactivator 1 (Mucoepidermoid carcinoma translocated protein 1) (Transducer of regulated cAMP response element-binding protein 1) (TORC-1) (Transducer of CREB protein 1), CRTC1, KIAA0616 MECT1 TORC1 WAMTP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CRTC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.