



TRNP1 Antibody, Biotin conjugated

Product Code	CSB-PA761223LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NT89
Immunogen	Recombinant Human TMF-regulated nuclear protein 1 protein (10-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TMF-regulated nuclear protein 1, TRNP1, C1orf225 TRNP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRNP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.