



SLC13A5 Antibody, Biotin conjugated

Product Code	CSB-PA768239LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YT5
Immunogen	Recombinant Human Solute carrier family 13 member 5 protein (156-204AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 13 member 5 (Na(+)/citrate cotransporter) (NaCT) (Sodium-coupled citrate transporter) (Sodium-dependent citrate transporter), SLC13A5, NACT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC13A5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.