




gelE Antibody

Product Code	CSB-PA769310LA01ELW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q833V7
Immunogen	Recombinant Enterococcus faecalis Gelatinase protein (193-510AA)
Raised In	Rabbit
Species Reactivity	Enterococcus faecalis
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterococcus faecalis
Research Area	Others
Target Names	gelE

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: gelE antibody at 2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 47.5 kDa

observed band size: 53 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.