



gelE Antibody, HRP conjugated

Product Code	CSB-PA769310LB01ELW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q833V7
Immunogen	Recombinant Enterococcus faecalis Gelatinase protein (193-510AA)
Raised In	Rabbit
Species Reactivity	Enterococcus faecalis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gelatinase (EC 3.4.24.30) (Coccolysin), gelE
Immunogen Species	Enterococcus faecalis (strain ATCC 700802 / V583)
Research Area	Others
Target Names	gelE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.